


# PA20 collimator: design of a multi-purpose collimator for SANS and GISANS at LLB

**Sylvain Désert**

Instrument Development Group  
CEA Saclay, France



40 m "Pinhole" geometry SANS / GISANS

Velocity selector 3,5-20 Å

Polarized neutron ( $> 4$  Å)

16 m collimator

19 m detector tank

2 sets of  $^3\text{He}$  detectors


Wide Q range ( $5 \cdot 10^{-4} - 0,5 \text{ \AA}^{-1}$ )

# Collimator

16 m long  
4 sections  
2 supports  
Vacuum

Rectangular apertures  
6 collimation lengths

Attenuator wheel  
Telescopic nose  
Inspection hatches



Foreseen “standard” collimation :  
 $25 \times 25 \text{ mm}^2$  entrance and  $12.5 \times 12.5 \text{ mm}^2$  exit  
i.e. full use of beam width

**For surface studies**


Slit shape beam

Hor. Or Ver. incident angle on the sample


Same options as for SANS mode

X


## Reflectometry on nanostructured $\text{CoSiO}_2$ / Si wafer




presence of weak stripe domains in the superferromagnetic phase



$\lambda = 7 \text{ \AA}$   
 Collimation : 2m  
 Diaphragmes : 1mm\*25mm



$Q_{\min}$ (1/Å)	$\lambda$ (Å)	Gain
$4 \cdot 10^{-4}$	9,2	13
$6 \cdot 10^{-4}$	9,2	8
$8 \cdot 10^{-4}$	6,5	5
$10^{-3}$	8,4	3



+ Dedicated free path (multibeam, Söller, reflective, ...)


## For studies of magnetic properties of samples

Polarization must be guided to the sample for 4+ Å neutrons


Constraints:

- Fit the guides and diaphragms
- No magnetic parts
- Telescopic part

Soft iron plates




Magnets  
(NdFeB, 0.5 T)




B (Gauss)

40  
39  
38  
37  
36  
35  
34  
33  
32

# Guide field







⇒ 35 G vertical ( $B_y$ ) and 2 G (along 30 mm) longitudinal ( $B_z$ )  
Adiabatic condition fulfilled


# Guide translations

85(H)x25(W) mm<sup>2</sup> inner section


3 elements of 4 m, and 1 element of 2 m


Translation length = 200 mm

⇒ for Guide IN, Guide OUT or free space IN





# Guide translations









# Blinds / Absorber




# Telescopic nose & Guide field



# Collimator support





# Collimator support





# 19 m collimation

Collimation length = Sample to detector





# Acknowledgements



## LLB Staff

Design office

P. Permingeat

P. Lavie

Polarization

S. Klimko

Motion control

F. Coneggio

P. Lambert

W. Josse

G. Koskas

Technicians

M. Detrez


A. Helary

Scientists


G. Chaboussant

A. Brûlet









Spacing : 100 nm  
Pore size : 40 nm



## GISANS on nanostructured $\text{CoSiO}_2$ / Si wafer



« domain » formation  
length scale : 125 nm





## Reflectometry on a Si wafer


$\lambda = 7 \text{ \AA}$

Collimation : 2m

Diaphragmes : 1mm\*25mm



## Reflectometry on nanostructured $\text{CoSiO}_2$ / Si wafer




presence of weak stripe domains in the superferromagnetic phase

$\lambda = 7 \text{ \AA}$

Collimation : 2m



Diaphragmes : 1mm\*25mm

PA20 - may 2015 - Co-SiO<sub>2</sub> granular / Si



Only 2 supports located at 25% of the length to minimize bending

(2 temporary supports needed during assembly as collimator is made out of 4 x 7.5 m pieces)





Instrument PA20

Collimator

Scientific modes:


- SANS
- GISANS
- VSANS

Polarization

Technical aspects








Rectangular apertures with slits

6 collimation lengths : 19, 16, 12, 8, 4, 2 m

Foreseen “standard” collimation :  $25 \times 25\text{ mm}^2$  entrance and  $12.5 \times 12.5\text{ mm}^2$  exit  
i.e. full use of beam width

Adiabatic rotation without loss of polarization:



$$\omega_L > 22. \omega_B \quad @ \lambda = 4 \text{ \AA}$$